



# Sterimatic Worldwide Ltd.

Griffin Mill, London Road, Stroud, Glos, GL5 2AZ, UK  
Tel: +44(0) 1453884944 Fax: +44(0) 1453886481 Mobile +44 (0) 7973361593  
Vat Reg: 800 8086 54

Registered Office: Greenmarley, Pinfarthings, Stroud, Glos, GL5 5JQ. Reg. No.: 4592319  
E-mail: [prowley@sterimatic.com](mailto:prowley@sterimatic.com) Web-site: [www.sterimatic.com](http://www.sterimatic.com)

Sterimatic Worldwide Ltd.

Sterimatic does not manufacture syringes, however there are many compatible types.

## Plastic multi-dose Syringes that are compatible with the Sterimatic system.

Simply remove needle locking nut from syringe and replace with Sterimatic sleeve.

### Simcrotech;

Standard type with draw-off-tube  
(1,2,5 or 10ml)



Bottle mounted, (2 or 5ml)



### Instrument supplies

Standard type with draw-off-tube  
(1,2,5 or 10ml)



Bottle mounted, (2 or 5ml)



As well as standard syringes that can be used by removing locking nut. Instrument supplies also supply a syringe with a metal luer-lock especially for use with Sterimatic.

### Phillips

Standard type(1,2,5 or 10ml)  
Or Bottle mount (2 or 5ml)



### Henke- Sass Ecomatic

Standard type or bottle mount (2 & 5ml)  
**NOTE: Henke - Sterimatic compatible  
Syringes are available only on request**



Phillips push-syringe.





Sterimatic Worldwide Ltd.

# Sterimatic Worldwide Ltd.

Griffin Mill, London Road, Stroud, Glos, GL5 2AZ, UK

Tel: +44(0) 1453884944 Fax: +44(0) 1453886481 Mobile +44 (0) 7973361593

Vat Reg: 800 8086 54

Registered Office: Greenmarley, Pinfarthings, Stroud, Glos, GL5 5JQ. Reg. No.: 4592319

E-mail: [prowley@sterimatic.com](mailto:prowley@sterimatic.com) Web-site: [www.sterimatic.com](http://www.sterimatic.com)

## TU – Trade-Unlimited

Standard type or bottle mount (2 & 5ml)

Or Bottle mount (2 or 5ml)

**NOTE: TU - Sterimatic compatible**

**Syringes are available only on request**



## Primatech

Standard type, (1,2,5 or 10ml)

Or Bottle mount (2 or 5ml)



## **Metal syringes**

The Sterimatic system will fit onto metal syringes such as the Kaycee, Socorex, Hauptner, Phillips & Henke-Sass. To enable this an adapter is usually required.

The adapters are available from the syringe manufacturer or from Sterimatic.

## **Phillips code: Sterimatic adapter**

(Sterimatic: ANA002)



## **Repeater syringe adapter: adapter SN402S**

(Sterimatic: ANA018)



## **Socorex code: 186.562**

(Sterimatic: ANA008)



## **Henke-Sass Roux revolver**

## **Henke code: Sterimatic adapter**

(Sterimatic: ANA014)



## **Hauptner codes: 16366.100 & 16102.250 & 16360200**

(Sterimatic: ANA015 & 16 & 17)



## **Kaycee code: 10145 + ANA002**

(Sterimatic: ANA007 & ANA002)





Sterimatic Worldwide Ltd.

# **Sterimatic Worldwide Ltd.**

Griffin Mill, London Road, Stroud, Glos, GL5 2AZ, UK

Tel: +44(0) 1453884944 Fax: +44(0) 1453886481 Mobile +44 (0) 7973361593

Vat Reg: 800 8086 54

Registered Office: Greenmarley, Pinfarthings, Stroud, Glos, GL5 5JQ. Reg. No.: 4592319

E-mail: [prowley@sterimatic.com](mailto:prowley@sterimatic.com) Web-site: [www.sterimatic.com](http://www.sterimatic.com)

Syringe type	Manufacturers code	Sterimatic code	Ex works Price from Sterimatic
Phillips 2&5ml tube fed Phillips repeater	Sterimatic adapter Sterimatic adapter	ANA002 ANA018	GBP 3.00 GBP 2.20
Socorex	186.562	ANA008	GBP 10.00
Henke-Sass Roux revolver	Sterimatic adapter	ANA014	GBP 10.20
Hauptner revolver	16366.100 & 16102.250 & 16360200	ANA015 & 16 & 17	Total GBP 10.20
Kaycee	10145 + ANA002	ANA007 & ANA002	Total GBP 9.85

Sterimatic does not add any margin to these items, however with delivery and currency changes costs it would be slightly cheaper to buy direct from the manufacturers.